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BIOLOG

Gram Positive Identification Test Panel

GP2 MicroPlate™

A1 Water	A2 α-Cyclodextrin	A3 β-Cyclodextrin	A4 Dextrin	A5 Glycogen	A6 Inulin	A7 Mannan	A8 Tween 40	A9 Tween 80	A10 N-Acetyl-D-Glucosamine	A11 N-Acetyl-β-D-Mannosamine	A12 Amygdalin
B1 L-Arabinose	B2 D-Arabitol	B3 Arbutin	B4 D-Cellobiose	B5 D-Fructose	B6 L-Fucose	B7 D-Galactose	B8 D-Galacturonic Acid	B9 Gentiobiose	B10 D-Gluconic Acid	B11 α-D-Glucose	B12 m-Inositol
C1 α-D-Lactose	C2 Lactulose	C3 Maltose	C4 Maltotriose	C5 D-Mannitol	C6 D-Mannose	C7 D-Melezitose	C8 D-Melibiose	C9 α-Methyl-D-Galactoside	C10 β-Methyl-D-Galactoside	C11 β-Methyl Glucose	C12 α-Methyl-D-Glucoside
D1 β-Methyl-D-Glucoside	D2 α-Methyl-D-Mannoside	D3 Palatinose	D4 D-Paicoose	D5 D-Raffinose	D6 L-Rhamnose	D7 D-Ribose	D8 Salicin	D9 Sedoheptulosan	D10 D-Sorbitol	D11 Stachyose	D12 Sucrose
E1 D-Tagatose	E2 D-Trehalose	E3 Turranose	E4 Xylitol	E5 D-Xylose	E6 Acetic Acid	E7 α-Hydroxybutyric Acid	E8 β-Hydroxybutyric Acid	E9 γ-Hydroxybutyric Acid	E10 p-Hydroxy-Phenylacetic Acid	E11 α-Ketoglutaric Acid	E12 α-Ketovaleric Acid
F1 Lactamide	F2 D-Lactic Acid Methyl Ester	F3 L-Lactic Acid	F4 D-Malic Acid	F5 L-Malic Acid	F6 Pyruvic Acid Methyl Ester	F7 Succinic Acid Mono-methyl Ester	F8 Propionic Acid	F9 Pyruvic Acid	F10 Succinamic Acid	F11 Succinic Acid	F12 N-Acetyl-L-Glutamic Acid
G1 L-Alaninamide	G2 D-Alanine	G3 L-Alanine	G4 L-Alanyl-Glycine	G5 L-Asparagine	G6 L-Glutamic Acid	G7 Glycyl-L-Glutamic Acid	G8 L-Proglutamic Acid	G9 L-Serine	G10 Putrescine	G11 2,3-Butanediol	G12 Glycerol
H1 Adenosine	H2 2'-Deoxy Adenosine	H3 Inosine	H4 Thymidine	H5 Uridine	H6 Adenosine-5'-Monophosphate	H7 Thymidine-5'-Monophosphate	H8 Uridine-5'-Monophosphate	H9 D-Fructose-6-Phosphate	H10 α-D-Glucose-1-Phosphate	H11 D-Glucose-6-Phosphate	H12 D,L-Glycerol Phosphate

FIGURE 1. Carbon Sources in GP2 MicroPlate

INTRODUCTION

The Biolog GP2 MicroPlate (Figure 1) is designed for identification and characterization of a very wide range of aerobic gram-positive bacteria. Biolog's MicroPlates and databases were first introduced in 1989, employing a novel, patented redox chemistry. This chemistry, based on reduction of tetrazolium, responds to the process of metabolism (i.e. respiration) rather than to metabolic by-products (e.g. acid). Biolog's chemistry works as a universal reporter of metabolism and simplifies the testing process as color developing chemicals do not need to be added. Since the GP2 MicroPlate is not dependent upon growth to produce identifications, it provides superior capability for all types of gram positive organisms: cocci, rods, and spore-forming rods all are identified with a single panel. The database for the GP2 MicroPlate is now over 310 species. It is by far the largest kit-based identification database available.

GP2 MICROPLATE

The Biolog GP2 MicroPlate performs 95 discrete tests simultaneously and gives a characteristic reaction pattern called a "metabolic fingerprint". These fingerprint reaction patterns provide a vast amount of information conveniently contained in a single Biolog MicroPlate. The metabolic fingerprint patterns are compared and identified using the MicroLog™ database software.

Other aerobic kit-based identification methods rely on a much smaller number of tests. Consequently, the significant limitation of these products is the limited number of species and organism types that they can identify. Furthermore, these products were designed to address the needs of routine clinical/hospital testing. The Biolog GP2 MicroPlate was designed to address the needs of a much wider range of users including environmental testing labs and animal and plant disease labs as well as clinical reference labs.

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Key : <X>: positive; <X-: mismatched positive; X: negative; X+: mismatched negative
 {X}: borderline; -X: less than A1 well

Color	1	2	3	4	5	6	7	8	9	10	11	12
A	0	47	< 114>< 163>	54	19	31	< 260>< 333>< 174>	47	48			
B	{ 83}	-1	{ 62} < 174>< 211>	19	{ 61}	41	< 132>{ 82} < 264>	53				
C	4	2	{ 92} 44 < 144>< 194>	50	{ 72}	36	18 < 147>< 112>					
D	39	-9	< 109>< 119>{ 70}	35	< 145>{ 75} { 71}	< 139>	37	< 198>				
E	15	< 183>{ 85}	23 < 146>	25	-11	36	46	15	47	{ 58}		
F	1	26	37 -24 { 69} { 71} { 56}	-7	< 138>	38	{ 58}	44				
G	16	23	{ 61} { 57} { 83} { 77}	21	{ 57} { 67}	41	{ 82}	< 106>				
H	< 117>< 103>< 112>< 122>< 134>	13	12	17	16	8	26	{ 78}				

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BIOLOG**Gram Negative Identification Test Panel****GN2 MicroPlate™**

A1 Water	A2 α-Cyclodextrin	A3 Dextrin	A4 Glycogen	A5 Tween 40	A6 Tween 80	A7 N-Acetyl-D-Galactosamine	A8 N-Acetyl-D-Glucosamine	A9 Adonitol	A10 L-Arabinose	A11 D-Arabitol	A12 D-Cellobiose
B1 L-Erythritol	B2 D-Fructose	B3 L-Fucose	B4 D-Galactose	B5 Gentiobiose	B6 α-D-Glucose	B7 m-Inositol	B8 α-D-Lactose	B9 Lactulose	B10 Maltose	B11 D-Mannitol	B12 D-Mannose
C1 D-Melibiose	C2 β-Methyl-D-Glucoside	C3 D- Psicose	C4 D-Raffinose	C5 L-Rhamnose	C6 D-Sorbitol	C7 Sucrose	C8 D-Trehalose	C9 Turannose	C10 Xylitol	C11 Pyruvic Acid Methyl Ester	C12 Succinic Acid Mono-Methyl-Ester
D1 Acetic Acid	D2 Cis-Aconitic Acid	D3 Citric Acid	D4 Formic Acid	D5 D-Galactonic Acid Lactone	D6 D-Galacturonic Acid	D7 D-Gluconic Acid	D8 D-Glucosaminic Acid	D9 D-Glucuronic Acid	D10 α-Hydroxybutyric Acid	D11 β-Hydroxybutyric Acid	D12 γ-Hydroxybutyric Acid
E1 p-Hydroxy Phenylacetic Acid	E2 Itaconic Acid	E3 α-Keto Butyric Acid	E4 α-Keto Glutaric Acid	E5 α-Keto Valeric Acid	E6 D,L-Lactic Acid	E7 Malonic Acid	E8 Propionic Acid	E9 Quinic Acid	E10 D-Saccharic Acid	E11 Sebacic Acid	E12 Succinic Acid
F1 Bromosuccinic Acid	F2 Succinamic Acid	F3 Glucuronamide	F4 L-Alaninamide	F5 D-Alanine	F6 L-Alanine	F7 L-Alanyl-glycine	F8 L-Asparagine	F9 L-Aspartic Acid	F10 L-Glutamic Acid	F11 Glycyl-L-Aspartic Acid	F12 Glycyl-L-Glutamic Acid
G1 L-Histidine	G2 Hydroxy-L-Proline	G3 L-Leucine	G4 L-Ornithine	G5 L-Phenylalanine	G6 L-Proline	G7 L-Pyroglutamic Acid	G8 D-Serine	G9 L-Serine	G10 L-Threonine	G11 D,L-Carnitine	G12 γ-Amino Butyric Acid
H1 Urocanic Acid	H2 Inosine	H3 Uridine	H4 Thymidine	H5 Phenethyl-amine	H6 Putrescine	H7 2-Aminoethanol	H8 2,3-Butanediol	H9 Glycerol	H10 D,L-α-Glycerol Phosphate	H11 α-D-Glucose-1-Phosphate	H12 D-Glucose-6-Phosphate

FIGURE 1. Carbon Sources in GN2 MicroPlate

INTRODUCTION

The Biolog GN2 MicroPlate (Figure 1) is designed for identification and characterization of a very wide range of aerobic gram-negative bacteria. Biolog's MicroPlates and databases were first introduced in 1989, employing a novel, patented redox chemistry. This chemistry, based on reduction of tetrazolium, responds to the process of metabolism (i.e. respiration) rather than to metabolic by-products (e.g. acid). Biolog's chemistry works as a universal reporter of metabolism and simplifies the testing process as color developing chemicals do not need to be added. Since the GN2 MicroPlate is not dependent upon growth to produce identifications, it provides superior capability for all types of gram negative organisms: fermenters, non-fermenters, and fastidious organisms all are identified with a single panel. The database for the GN2 MicroPlate is now over 500 species. It is by far the largest kit-based identification database available.

GN2 MICROPLATE

The Biolog GN2 MicroPlate performs 95 discrete tests simultaneously and gives a characteristic reaction pattern called a "metabolic fingerprint". These fingerprint reaction patterns provide a vast amount of information conveniently contained in a single Biolog MicroPlate. The metabolic fingerprint patterns are compared and identified using the MicroLog™ database software.

Other aerobic kit-based identification methods rely on a much smaller number of tests. Consequently, the significant limitation of these products is the limited number of species and organism types that they can identify. Furthermore, these products were designed to address the needs of routine clinical/hospital testing. The Biolog GN2 MicroPlate was designed to address the needs of a much wider range of users including environmental testing labs and animal and plant disease labs as well as clinical reference labs.

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Key : <X>: positive; <X-: mismatched positive; X: negative; X+: mismatched negative
{X}: borderline; -X: less than A1 well

Color	1	2	3	4	5	6	7	8	9	10	11	12
A	0	-7	{ 65}	< 111>< 107>{ 71}	-7	12	-3	< 271>	31	56		
B	-4	< 107>	-13	27	9	< 250>	-12	-11	-1	-9	< 159>{ 79}	
C	30	-8	49	-5	-12	3	-4	30	30	-5	< 100>< 104>	
D	{ 78}	< 230>< 342>< 134>	-15	-1	< 256>	-11	2	{ 71}	< 162>	-1		
E	< 275>< 207>{ 83}	< 261>< 153>< 170>	-9	< 203>< 146>	2	-2	< 162>					
F	< 277>	36	13	{ 81}	< 181>< 161>< 122>< 464>< 355>< 307>	47	57					
G	3	< 109>	25	< 142>	-3	< 158>< 228>	24	< 185>< 121>{ 80}	< 194>			
H	12	< 131>	15	5	1	< 142>{ 80}	27	< 141>	36	0	{ 67}	

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
Identification flow charts

Differentiation via Gram stains and cell morphology

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Copyright© by Chiang Mai University
All rights reserved